

Declaration of Conformity

Medema-A/S declares on our sole responsibility that the product mentioned below are manufactured in accordance with our incorporated standards and are in conformity with the following directives and standards if it is used according to our instructions:

1506-2051 BOBIKE ONE BRACKET KIT

Can be used to mount on a Mini Crosser and function as a luggage carrier as if it was mounted on a bicycle under the following circumstances.

The scooter must be fitted with anti-tip devices.

The seat of the scooter must always be in the lowest position.

The maximum user weight of 175 kg is including driver and child all together.

The child alone may not way more than 22 kg as specified in the user manual for BoBike one child seat.

The program settings of the scooter controller must be ensured are changed to the following values:

Forward acceleration 40

Reverse deceleration 25

The installation must be performed by Mini Crosser dealer or technician.

Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

EN 12184:2014 Electrically powered wheelchairs, scooters and their chargers - Requirements and test methods.

EN 11243:2016 Cycles – Luggage carriers for bicycles – Requirements and test methods.

Directive 2006/42/EC of the European Parliament and of the Council of 17 May 2006 on machinery, and amending Directive 95/16/EC.

ISO 14971:2012 Medical devices – Application of risk management to medical devices.

Title: General Manager

Name: Klaus Smedegaard



Date: 04-10-2021 Signature

Medema A/S

Enggårdvej 7
Snebjerg
7400 Herning

CVR/VAT: DK 51 89 27 12
SRN: DK-MF-000000981

Phone: 70 10 17 55

info@medema.com
www.medema.com



medemagruppen