

# Declaration of Conformity

We, the manufacturer:

Medema A/S  
Enggårdvej 7  
DK-7400 Herning

CVR/VAT: DK 51 89 27 12

Hereby declares that the product:

## **1506-2051 BOBIKE ONE BRACKET KIT**

Complies with the Medical Device Directive 93/42/EEC when mounted on the Mini Crosser scooter M-model and X-model under the following circumstances:

The scooter must be fitted with anti-tip devices.

The seat of the scooter must always be in the lowest position.

The maximum user weight of 175 kg is including driver and child all together.

The child alone may not weigh more than 22 kg as specified in the user manual for BoBike one child seat.

The program settings of the scooter controller must be ensured are changed to the following values:

Forward acceleration 40  
Reverse deceleration 25

The installation must be performed by Mini Crosser dealer or technician.

Title: R&D Manager

Name: Klaus Smedegaard

Date: 21-03-2019

Signature: \_\_\_\_\_



Medema A/S

Enggårdvej 7  
7400 Herning

T +45 70 10 20 54

F + 45 97 16 85 82

[www.medema.com](http://www.medema.com)



medemagruppen